

Product Datasheet

Human sRANKL Protein (His Tag) (orb80331)

Catalog Number	orb80331
Description	Recombinant of Human sRANKL protein (His Tag)
Preservatives	sRANKL protein solution (0.5mg/ml) containing 20mM Tris-HCl buffer (pH8.0), 20% glycerol, 0.1M NaCl and 1mM DTT.
Form/Appearance	Sterile Filtered colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Cytokines And Growth Factors
Protein Sequence	MGSSHHHHHH SSGLVPRGSH MIRA EKAMVD GSWLDLAKRS KLEAQPF AHL TINATDIPSG SHKVSLS SWY HDRGWAKISN MTF SNGKLIV NQDGFYYLYA NICFRHHETS GD LATEYLQL M VYVTKTSIK IPSSHTLMKG GSTKYWSGNS EFHFYSINVG GFFKLRS GEE ISIEVSNPSL LDPDQDATYF GAFKVRDID
Purity	Greater than 80.0% as determined by SDS-PAGE.
Source	Escherichia Coli
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)