

Product Datasheet

Streptavidin-NC Protein (orb81397)

Catalog Number	orb81397
Description	Recombinant of Streptavidin-NC protein
Preservatives	The sterile solution contains 10mM K ₂ HPO ₄ -KH ₂ PO ₄ , pH 7.3.
Form/Appearance	Sterile Liquid formulation.
Storage	Stability: Streptavidin-NC although stable at 4°C for 3 weeks, should be stored below -18°C. Please prevent freeze thaw cycles
Note	For research use only
Application notes	Recombinant & Natural Proteins
Purity	Greater than 93.0% as determined by SDS-PAGE.
Source	Escherichia Coli
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)