

Product Datasheet

ubiG Antibody (orb834244)

Catalog Number	orb834244
Category	Antibodies
Description	This ubiG Antibody is an unconjugated polyclonal product. It targets ubiG using a synthesized peptide derived from Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC) ubiG as the immunogen. This antibody is suitable for DOT, ELISA. Purification: Antigen Affinity.
Target	ubiG
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	Unconjugated
Reactivity	E. coli
Form/Appearance	Liquid
Buffer/Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Purification	Antigen Affinity
Immunogen	Recombinant Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC) ubiG protein
UniProt ID	Q8FFP0
Tested applications	ELISA, WB

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Application notes

1. Biorbyt guarantees the antibody purity > 95% confirmed by SDS-PAGE. 2. Biorbyt guarantees the antibody titer > 1: 64,000 confirmed by ELISA. 3. Biorbyt guarantees to offer antigen-specific Western Blot validation data. 4. Biorbyt guarantees 100% risk-free: No qualified antibody, no charge.

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)