

Product Datasheet

Human C3a anaphylatoxin ELISA Kit (orb865345)

Catalog Number	orb865345
Category	Assays and Kits
Description	Human C3a anaphylatoxin ELISA Kit
Target	Complement C3
Reactivity	Human
Immunogen	Expression system for standard: E.coli; Immunogen sequence: S672-R748
UniProt ID	P01024
Tested applications	ELISA
Specificity	Natural and recombinant human C3
Cross Reactivity	There is no detectable cross-reactivity with other relevant proteins.
Range	31.2 pg/ml - 2,000 pg/ml
Sensitivity	20 pg/ml
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA, citrate) and urine
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.


7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



Human C3a anaphylatoxin PicoKine ELISA Kit Standard Curve.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)