

## Product Datasheet

### EB3 Rabbit Monoclonal Antibody (orb865900)

|                     |                                                                                                                                   |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------|
| Catalog Number      | orb865900                                                                                                                         |
| Category            | Antibodies                                                                                                                        |
| Description         | EB3 Rabbit Monoclonal Antibody                                                                                                    |
| Clonality           | Monoclonal                                                                                                                        |
| Species/Host        | Rabbit                                                                                                                            |
| Isotype             | IgG                                                                                                                               |
| Conjugation         | Unconjugated                                                                                                                      |
| Reactivity          | Human, Mouse, Rat                                                                                                                 |
| Form/Appearance     | Liquid                                                                                                                            |
| Concentration       | 0.5mg/ml                                                                                                                          |
| Purification        | Affinity-chromatography                                                                                                           |
| Immunogen           | A synthesized peptide derived from human EB3                                                                                      |
| UniProt ID          | <b>Q9UPY8</b>                                                                                                                     |
| MW                  | 30 kDa                                                                                                                            |
| Tested applications | FC, ICC, IF, IHC, WB                                                                                                              |
| Application notes   | WB 1:500-2000IHC 1:50-200ICC/IF 1:50-200FC 1:50. Restore with deionized water (or equivalent) for reconstitution volume of 1.0 mL |
| Antibody Type       | Primary Antibody                                                                                                                  |

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

#### Biorbyt LLC


68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

**Clone Number** 18M56

**Storage** Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

**Note** For research use only



Western blot analysis of EB3 expression in Ratmuscle cell lysate.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558