

## Product Datasheet

### Anti-MED4 Rabbit Monoclonal Antibody (orb866729)

|                            |                                                                                                                                         |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>      | orb866729                                                                                                                               |
| <b>Description</b>         | Anti-MED4 Rabbit Monoclonal Antibody                                                                                                    |
| <b>Species/Host</b>        | Rabbit                                                                                                                                  |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                                       |
| <b>Conjugation</b>         | Unconjugated                                                                                                                            |
| <b>Tested Applications</b> | ICC, IF, IHC, WB                                                                                                                        |
| <b>Immunogen</b>           | A synthesized peptide derived from human MED4                                                                                           |
| <b>Form/Appearance</b>     | Liquid                                                                                                                                  |
| <b>Concentration</b>       | Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.                                            |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                                                                   |
| <b>Application notes</b>   | WB 1:500-1:2000IHC 1:50-1:200ICC/IF 1:50-1:200. Restore with deionized water (or equivalent) for reconstitution volume of 1.0 mL        |
| <b>Isotype</b>             | IgG                                                                                                                                     |
| <b>Clonality</b>           | Monoclonal                                                                                                                              |
| <b>Clone Number</b>        | 26M85                                                                                                                                   |
| <b>Antibody Type</b>       | Primary Antibody                                                                                                                        |


#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

|                        |                                 |
|------------------------|---------------------------------|
| <b>MW</b>              | 35 kDa                          |
| <b>Uniprot ID</b>      | <b>Q9NPJ6</b>                   |
| <b>Expiration Date</b> | 12 months from date of receipt. |



Western blot analysis of MED4 expression in HepG2 cell lysate.

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)