

Product Datasheet

MOUSE ANTI HUMAN A2M:Biotin (orb867195)

Catalog Number	orb867195
Description	Mouse monoclonal antibody against A2M conjugated to Biotin
Species/Host	Mouse
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	E. coli derived recombinant human A2M (Ser24-Ala1474)
Target	A2M
Preservatives	Purified IgG conjugated to Biotin - liquid
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	IgG concentration 0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Clone Number	M5
Uniprot ID	P01023
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)