

Product Datasheet

Cynomolgus ROR2 protein (orb867373)

Catalog Number	orb867373
Category	Proteins
Description	Recombinant Cynomolgus ROR2 (C-hFc) is produced by Human Cells. The target gene encoding Val35-Gly403 is expressed with a C-hFc tag.
Tag	C-Human Fc
UniProt ID	A0A2K5UT30
Expression System	Human Cells
Endotoxins	Less than 1 EU/ug
Biological Origin	Cynomolgus
Expression Region	Val35-Gly403
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)