

## Product Datasheet

### Mouse Anti-Rabbit IgG H&L, Cy5.5 conjugated (orb868500)

|                             |                                                                                                                                                              |
|-----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>       | orb868500                                                                                                                                                    |
| <b>Category</b>             | Antibodies                                                                                                                                                   |
| <b>Description</b>          | The Mouse Anti-Rabbit IgG H&L, Cy5.5 conjugated is a Cy5.5 conjugated polyclonal antibody suitable for FC, ICC, IF. This product reacts with Rabbit samples. |
| <b>Clonality</b>            | Polyclonal                                                                                                                                                   |
| <b>Species/Host</b>         | Mouse                                                                                                                                                        |
| <b>Isotype</b>              | IgG                                                                                                                                                          |
| <b>Conjugation</b>          | Cy5.5                                                                                                                                                        |
| <b>Reactivity</b>           | Rabbit                                                                                                                                                       |
| <b>Form/Appearance</b>      | Liquid                                                                                                                                                       |
| <b>Concentration</b>        | 2.0 mg/ml                                                                                                                                                    |
| <b>Buffer/Preservatives</b> | 10 mM TBS (pH=7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% glycerol.                                                                                      |
| <b>Purification</b>         | Affinity purified by Protein G                                                                                                                               |
| <b>Immunogen</b>            | Native rabbit IgG                                                                                                                                            |
| <b>Tested applications</b>  | FC, ICC, IF                                                                                                                                                  |
| <b>Antibody Type</b>        | Secondary Antibody                                                                                                                                           |
| <b>Storage</b>              | Store at -20°C for one year. Avoid repeated freeze/thaw cycles.                                                                                              |

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

|                        |                                 |
|------------------------|---------------------------------|
| <b>Note</b>            | For research use only           |
| <b>Expiration Date</b> | 12 months from date of receipt. |

---

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)