

Product Datasheet

Rabbit Anti-Rat IgG H&L, Cy5 conjugated (orb868522)

Catalog Number	orb868522
Description	Rabbit Anti-Rat IgG H&L, Cy5 conjugated
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Cy5
Tested Applications	FC, IF
Immunogen	Native Rat IgG
Preservatives	10 mM TBS (pH=7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% glycerol.
Concentration	2.0 mg/ml
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Antibody Type	Secondary Antibody
Dilution Range	IF=1:100-1000, Flow-Cyt=1:100-1000
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)