

Product Datasheet

Goat Anti-Rabbit IgG H&L, Gold conjugated (orb868897)

Catalog Number	orb868897
Category	Antibodies
Description	Goat Anti-Rabbit IgG H&L, Gold conjugated
Clonality	Polyclonal
Species/Host	Goat
Isotype	IgG
Conjugation	Gold
Reactivity	Rabbit
Concentration	0.4 mg/ml
Buffer/Preservatives	20mM TBS (pH8.0) with 1% rAlbumin and 0.02% Proclin300.
Purification	Affinity purified by Protein G, nonspecific adsorbed
Immunogen	Native rabbit IgG
Tested applications	ICA, IEM
Dilution range	ICA=1:20-200, IEM=1:20-200
Antibody Type	Secondary Antibody
Storage	Store at 2-8°C for 3-6 months. Avoid repeated freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)