

## Product Datasheet

### Mouse Human IgG H&L Monoclonal Antibody (Biotin) (orb868930)

|                             |                                                                                                                                                                                                   |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>       | orb868930                                                                                                                                                                                         |
| <b>Category</b>             | Antibodies                                                                                                                                                                                        |
| <b>Description</b>          | The Mouse Anti-Human IgG H&L (B5F5) mAb, Biotin conjugated is a Biotin conjugated monoclonal antibody suitable for ELISA, FC, ICC, IF, IHC-Fr, IHC-P, WB. This product reacts with Human samples. |
| <b>Clonality</b>            | Monoclonal                                                                                                                                                                                        |
| <b>Species/Host</b>         | Mouse                                                                                                                                                                                             |
| <b>Isotype</b>              | IgG1                                                                                                                                                                                              |
| <b>Conjugation</b>          | Biotin                                                                                                                                                                                            |
| <b>Reactivity</b>           | Human                                                                                                                                                                                             |
| <b>Form/Appearance</b>      | Liquid                                                                                                                                                                                            |
| <b>Concentration</b>        | 2.0 mg/ml                                                                                                                                                                                         |
| <b>Buffer/Preservatives</b> | 10 mM TBS (pH=7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% glycerol.                                                                                                                           |
| <b>Purification</b>         | Affinity purified by Protein G                                                                                                                                                                    |
| <b>Immunogen</b>            | Native Human IgG                                                                                                                                                                                  |
| <b>Tested applications</b>  | ELISA, FC, ICC, IF, IHC-Fr, IHC-P, WB                                                                                                                                                             |
| <b>Antibody Type</b>        | Secondary Antibody                                                                                                                                                                                |
| <b>Clone Number</b>         | B5F5                                                                                                                                                                                              |

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

|                        |                                                                 |
|------------------------|-----------------------------------------------------------------|
| <b>Storage</b>         | Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| <b>Note</b>            | For research use only                                           |
| <b>Expiration Date</b> | 12 months from date of receipt.                                 |

---

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)