

Product Datasheet

Goat Anti-Rat IgG H&L, HRP conjugated (orb868939)

Catalog Number	orb868939
Category	Antibodies
Description	Goat Anti-Rat IgG H&L, HRP conjugated
Clonality	Polyclonal
Species/Host	Goat
Isotype	IgG
Conjugation	HRP
Reactivity	Rat
Concentration	2.0 mg/ml
Buffer/Preservatives	10 mM TBS (pH=7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% glycerol.
Purification	Affinity purified by Protein G
Immunogen	Native Rat IgG
Tested applications	ELISA, IHC-Fr, IHC-P, WB
Dilution range	WB=1:2000-20000, IHC-P=1:100-1000, IHC-F=1:100-1000, ELISA=1:2000-20000
Antibody Type	Secondary Antibody
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558