

## Product Datasheet

# Phospho-SGK1 (Thr256) Rabbit Polyclonal Antibody (PE-Cy5) (orb886552)

|                            |   |
|----------------------------|---|
| <b>Catalog Number</b>      | orb886552   |
| <b>Description</b>         | Phospho-SGK1 (Thr256) Rabbit Polyclonal Antibody (PE-Cy5)   |
| <b>Species/Host</b>        | Rabbit  |
| <b>Reactivity</b>          | Human, Mouse  |
| <b>Conjugation</b>         | PE/Cy5  |
| <b>Tested Applications</b> | FC, ICC, IF   |
| <b>Immunogen</b>           | KLH conjugated Synthesised phosphopeptide derived from human SGK1 around the phosphorylation site of Thr256 TS(p-T)FC |
| <b>Target</b>              | SGK1  |
| <b>Preservatives</b>       | 0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.  |
| <b>Concentration</b>       | 1mg/ml  |
| <b>Storage</b>             | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.                                       |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Antibody Type</b>       | Primary Antibody  |
| <b>Dilution Range</b>      | Flow-Cyt=2ug/Test, ICC/IF=1:50-200, IF=1:100-500  |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)