

Product Datasheet

Phospho-AMPK beta 1 + AMPK beta 2 (Ser182/Ser184) Rabbit Polyclonal Antibody (PE- Cy5) (orb889405)

Catalog Number	orb889405
Description	Phospho-AMPK beta 1 + AMPK beta 2 (Ser182/Ser184) Rabbit Polyclonal Antibody (PE-Cy5)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PE/Cy5
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated synthesised phosphopeptide derived from human AMPK beta 1 around the phosphorylation site of Ser182 SS(p-S)PP
Target	PRKAB1
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Antibody Type	Primary Antibody
Dilution Range	Flow-Cyt=2ug/Test, ICC/IF=1:50-200, IF=1:100-500

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Expiration Date

12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)