

## Product Datasheet

### Phospho-PAK6 (Ser560) Rabbit Polyclonal Antibody (PE-Cy7) (orb902630)

|                             |                                                                                                                                                                                                             |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>       | orb902630                                                                                                                                                                                                   |
| <b>Category</b>             | Antibodies                                                                                                                                                                                                  |
| <b>Description</b>          | Phospho-PAK6 (Ser560) Rabbit Polyclonal Antibody (PE-Cy7) is a PE/Cy7 conjugated antibody targeting PAK6. This antibody is suitable for FC, ICC, IF. It exhibits reactivity with Human, Mouse, Rat samples. |
| <b>Target</b>               | PAK6                                                                                                                                                                                                        |
| <b>Clonality</b>            | Polyclonal                                                                                                                                                                                                  |
| <b>Species/Host</b>         | Rabbit                                                                                                                                                                                                      |
| <b>Isotype</b>              | IgG                                                                                                                                                                                                         |
| <b>Conjugation</b>          | PE/Cy7                                                                                                                                                                                                      |
| <b>Reactivity</b>           | Human, Mouse, Rat                                                                                                                                                                                           |
| <b>Predicted Reactivity</b> | Bovine, Canine, Equine, Gallus, Porcine, Rabbit, Sheep                                                                                                                                                      |
| <b>Concentration</b>        | 1mg/ml                                                                                                                                                                                                      |
| <b>Buffer/Preservatives</b> | 0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.                                                                                                                                      |
| <b>Purification</b>         | Affinity purified by Protein A                                                                                                                                                                              |
| <b>Immunogen</b>            | KLH conjugated Synthesised phosphopeptide derived from human PAK6 around the phosphorylation site of Ser560 RK(p-S)LV                                                                                       |
| <b>Tested applications</b>  | FC, ICC, IF                                                                                                                                                                                                 |

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

|                       |                                                                                 |
|-----------------------|---------------------------------------------------------------------------------|
| <b>Dilution range</b> | ICC/IF=1:100-500, Flow-Cyt=1ug/Test, IF=1:100-500                               |
| <b>Antibody Type</b>  | Primary Antibody                                                                |
| <b>Storage</b>        | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| <b>Note</b>           | For research use only                                                           |

---

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)