

Product Datasheet

PCR-grade Water (orb93976)

Catalog Number	orb93976
Category	Tools
Description	PCR-grade Water is intended for use in molecular biology applications including PCR and RT-PCR. The ultra-pure and sterile filtered water is manufactured free of detectable inhibitors, contaminants or enzymatic activity.
Form/Appearance	liquid
Storage	store at 4 °C or -20 °C
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)