

Product Datasheet

T7 tag Mouse mAb, Cy5 conjugated (orb969297)

Catalog Number	orb969297
Category	Antibodies
Description	Mouse monoclonal antibody to T7-tag/Cy5
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG
Conjugation	Cy5
Reactivity	All
Concentration	1mg/ml
Buffer/Preservatives	Lyophilized or Liquid
Purification	Affinity purified by Protein G
Immunogen	KLH conjugated synthetic peptide T7-tag MASMTGGQQMG
Tested applications	FC, ICC, IF
Dilution range	IF=1:200-1000, Flow-Cyt=1ug/Test, ICC/IF=1:100-500
Antibody Type	Primary Antibody
Clone Number	3G7
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Expiration Date

12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)