

Product Datasheet

Mouse Anti-Rabbit IgG mAb (orb973625)

Catalog Number	orb973625
Category	Antibodies
Description	Mouse Anti-Rabbit IgG mAb
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG1
Conjugation	Unconjugated
Reactivity	Rabbit
Form/Appearance	Lyophilized
Buffer/Preservatives	0.01M PBS (pH7.4).
Purification	Affinity purified by Protein G
Immunogen	Native Rabbit IgG full length
Tested applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Dilution range	WB=1:5000-50000, IHC-P=1:100-500, IHC-F=1:1000-5000, ICC/IF=1:1000-5000, IF=1:1000-5000, ELISA=1:5000-50000
Antibody Type	Secondary Antibody
Clone Number	B6E33
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Note For research use only

Expiration Date 12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)