

## Product Datasheet

# Phospho-JNK1 + JNK2 + JNK3 (T183 + T183 + T221) Rabbit Polyclonal Antibody (Cy3) (orb975757)

|                             |   |
|-----------------------------|---|
| <b>Catalog Number</b>       | orb975757   |
| <b>Category</b>             | Antibodies  |
| <b>Description</b>          | Phospho-JNK1 + JNK2 + JNK3 (T183 + T183 + T221) Rabbit Polyclonal Antibody (Cy3) is a Cy3 conjugated antibody targeting MAPK8. This antibody is suitable for IF. It exhibits reactivity with Human, Mouse, Rat samples. |
| <b>Target</b>               | MAPK8   |
| <b>Clonality</b>            | Polyclonal  |
| <b>Species/Host</b>         | Rabbit  |
| <b>Isotype</b>              | IgG   |
| <b>Conjugation</b>          | Cy3   |
| <b>Reactivity</b>           | Human, Mouse, Rat   |
| <b>Predicted Reactivity</b> | Bovine, Canine, Porcine   |
| <b>Concentration</b>        | 1mg/ml  |
| <b>Buffer/Preservatives</b> | 0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.  |
| <b>Purification</b>         | Affinity purified by Protein A  |
| <b>Immunogen</b>            | KLH conjugated Synthesised phosphopeptide derived from human JNK1/JNK2/JNK3 around the phosphorylation site of T183/T183/T221 MM(p-T)PY   |

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

|                            |   |
|----------------------------|---|
| <b>Tested applications</b> | IF  |
| <b>Dilution range</b>      | IF=1:100-500  |
| <b>Antibody Type</b>       | Primary Antibody  |
| <b>Storage</b>             | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| <b>Note</b>                | For research use only   |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |

---

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)