

Product Datasheet

Collagen alpha-1 (I) chain Rabbit Polyclonal Antibody (Cy3) (orb977298)

Catalog Number	orb977298
Category	Antibodies
Description	Collagen alpha-1 (I) chain Rabbit Polyclonal Antibody (Cy3)
Target	COL1A1
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	Cy3
Reactivity	Rat
Predicted Reactivity	Human, Mouse
Concentration	1mg/ml
Buffer/Preservatives	Preservative: 0.02% Proclin300, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4.
Purification	affinity purified by Protein A
Immunogen	KLH conjugated synthetic peptide derived from human Collagen alpha-1(I) chain (1051-1218aa)
Tested applications	IF
Dilution range	IF=1:100-500

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Antibody Type	Primary Antibody
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)