

# Product Datasheet

## GABBR2 Rabbit Polyclonal Antibody (APC) (orb995206)

|                      |                                                                            |
|----------------------|----------------------------------------------------------------------------|
| Catalog Number       | orb995206                                                                  |
| Category             | Antibodies                                                                 |
| Description          | GABBR2 Rabbit Polyclonal Antibody (APC)                                    |
| Target               | GABBR2                                                                     |
| Clonality            | Polyclonal                                                                 |
| Species/Host         | Rabbit                                                                     |
| Isotype              | IgG                                                                        |
| Conjugation          | APC                                                                        |
| Reactivity           | Human, Mouse, Rat                                                          |
| Predicted Reactivity | Bovine, Canine, Equine, Gallus, Rabbit                                     |
| Concentration        | 1mg/ml                                                                     |
| Buffer/Preservatives | 0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.     |
| Purification         | Affinity purified by Protein A                                             |
| Immunogen            | KLH conjugated synthetic peptide derived from human GABBR2 (501-650/941aa) |
| Tested applications  | IF                                                                         |
| Dilution range       | IF=1:100-500                                                               |
| Antibody Type        | Primary Antibody                                                           |

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

**Storage**

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

**Note**

For research use only

---

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)