

Product Datasheet

Human PLA1A ELISA Kit orb405016

Description	Human Phospholipase A1 member A(PLA1A) ELISA kit
Reactivity	Human
Range	0.312 ng/mL-20 ng/mL
Sensitivity	0.078 ng/mL
Tested Applications	ELISA
Target	PLA1A
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q53H76
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates. • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm