

Product Datasheet

Human GPX4 ELISA Kit orb405019

Description	Human Phospholipid hydroperoxide glutathione peroxidase, mitochondrial(GPX4) ELISA kit
Reactivity	Human
Range	7.8 µU/mL-500 µU/mL
Sensitivity	1.95 µU/mL
Tested Applications	ELISA
Target	GPX4
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P36969
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm