

Product Datasheet

Human PLAC1 ELISA Kit (orb405025)

Description	Human Placenta-specific protein 1(PLAC1) ELISA kit
Reactivity	Human
Range	1.56 ng/mL-100 ng/mL
Sensitivity	0.39 ng/mL
Tested Applications	ELISA
Target	PLAC1
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9HBJ0
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm