
Product Datasheet

Human PZP ELISA Kit (orb405046)

Description	Human Pregnancy zone protein(PZP) ELISA kit
Reactivity	Human
Range	2 µg/mL-75 µg/mL
Sensitivity	2 µg/mL
Tested Applications	ELISA
Target	PZP
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P20742
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma and tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm