

Product Datasheet

Human CFP ELISA Kit orb405077

Description	Human Properdin(CFP) ELISA kit
Reactivity	Human
Range	12.5 ng/ml-800 ng/ml
Sensitivity	3.12ng/ml
Tested Applications	ELISA
Target	CFP
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, cell culture supernates, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm