

---

## Product Datasheet

### Human Interferon gamma ELISA Kit (orb405088)

<b>Description</b>	Human Proteasome activator complex subunit 1(PSME1) ELISA kit
<b>Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Target</b>	PSME1
<b>Note</b>	For research use only.
<b>Uniprot ID</b>	<b>Q06323</b>
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm