

Product Datasheet

Human FAM132B ELISA Kit orb405107

Description	Human Protein FAM132B(FAM132B) ELISA kit
Reactivity	Human
Range	15.6 pg/mL-1000 pg/mL
Sensitivity	3.9 pg/mL
Tested Applications	ELISA
Target	FAM132B
Note	For research use only.
Sample Types	serum, plasma, cell culture supernates, tissue homogenates, cell lysates
Uniprot ID	Q4G0M1
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, cell culture supernates, tissue homogenates, cell lysates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm