

Product Datasheet

Human FAM132B ELISA Kit orb405107

Description	Human Protein FAM132B(FAM132B) ELISA kit
Reactivity	Human
Range	15.6 pg/ml-1000 pg/ml
Sensitivity	3.9 pg/ml
Tested Applications	ELISA
Target	FAM132B
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, cell culture supernates, tissue homogenates, cell lysates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm