

Product Datasheet

Human FAM13A ELISA Kit orb405108

Description	Human Protein FAM13A(FAM13A) ELISA kit
Reactivity	Human
Range	23.44 pg/mL-1500 pg/mL
Sensitivity	5.86 pg/mL
Tested Applications	ELISA
Target	FAM13A
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	O94988
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm