

Product Datasheet

Human WNT2 ELISA Kit orb405146

Description	Human Protein Wnt-2(WNT2) ELISA kit
Reactivity	Human
Range	0.45 ng/ml-30 ng/ml
Sensitivity	0.11 ng/ml
Tested Applications	ELISA
Target	WNT2
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm