

Product Datasheet

Human WNT5B ELISA Kit orb405150

Description	Human Protein Wnt-5b(WNT5B) ELISA kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Sensitivity	0.039 ng/mL
Tested Applications	ELISA
Target	WNT5B
Note	For research use only.
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q9H1J7
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, cell culture supernates, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm