

Product Datasheet

Human ERBB3 ELISA Kit orb405210

Description	Human Receptor tyrosine-protein kinase erbB-3(ERBB3) ELISA kit
Reactivity	Human
Range	125 pg/mL-8000 pg/mL
Sensitivity	31.25 pg/mL
Tested Applications	ELISA
Target	ERBB3
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P21860
Expiration Date	Please enquire.
Application Notes	• Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm