

Product Datasheet

Human PTPRZ1 ELISA Kit orb405216

Description	Human Receptor-type tyrosine-protein phosphatase zeta(PTPRZ1) ELISA kit
Reactivity	Human
Range	93.75 pg/ml-6000 pg/ml
Sensitivity	23.44 pg/ml
Tested Applications	ELISA
Target	PTPRZ1
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm