

Product Datasheet

Human SFRS1 ELISA Kit orb405291

Description	Human Serine/arginine-rich splicing factor 1(SFRS1) ELISA kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Sensitivity	7.81 pg/mL
Tested Applications	ELISA
Target	SFRS1
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q07955
Expiration Date	Please enquire.
Application Notes	• Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm