

Product Datasheet

Human SAA2 ELISA Kit orb405314

Description	Human serum amyloid A2 (SAA2) ELISA kit
Reactivity	Human
Range	0.312 ng/ml-20 ng/ml
Sensitivity	0.078 ng/ml
Tested Applications	ELISA
Target	SAA2
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates.• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm