

Product Datasheet

Human SCUBE1 ELISA Kit orb405327

Description	Human Signal Peptide, CUB and EGF-like Domain-containing Protein 1(SCUBE1)ELISA Kit
Reactivity	Human
Range	0.625 ng/mL-40 ng/mL
Sensitivity	0.156 ng/mL
Tested Applications	ELISA
Target	SCUBE1
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q8IWIY4
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm