
Product Datasheet

Human SORT1 ELISA Kit (orb405350)

Description	Human Sortilin(SORT1) ELISA kit
Reactivity	Human
Range	46.88 pg/mL-3000 pg/mL
Sensitivity	11.72 pg/mL
Tested Applications	ELISA
Target	SORT1
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q99523
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates and cell lysates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm