

Product Datasheet

Human STATH ELISA Kit orb405364

Description	Human Statherin(STATH) ELISA kit
Reactivity	Human
Range	78 ng/mL-5000 ng/mL
Sensitivity	19.5 ng/mL
Tested Applications	ELISA
Target	STATH
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates, saliva
Uniprot ID	P02808
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates, saliva • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm