

Product Datasheet

Human **SULT1A1 ELISA Kit** orb405368

Description	Human Sulfotransferase 1A1(SULT1A1) ELISA kit
Reactivity	Human
Range	25 pg/mL-1600 pg/mL
Sensitivity	6.25 pg/mL
Tested Applications	ELISA
Target	SULT1A1
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P50225
Expiration Date	Please enquire.
Application Notes	• Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm