
Product Datasheet

Human CD247 ELISA Kit (orb405396)

Description	Human T-cell surface glycoprotein CD3 zeta chain (CD247) ELISA kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Sensitivity	0.039 ng/mL
Tested Applications	ELISA
Target	CD247
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P20963
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates, cell lysates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm