

## Product Datasheet

### Human TARS ELISA Kit orb405413

<b>Description</b>	Human Threonyl-tRNA synthetase, cytoplasmic(TARS) ELISA kit
<b>Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Target</b>	TARS
<b>Note</b>	For research use only.
<b>Uniprot ID</b>	<b>P26639</b>
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	• Sample Type: serum, plasma, cell lysates, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm