

Product Datasheet

Human TARS ELISA Kit orb405413

Description	Human Threonyl-tRNA synthetase, cytoplasmic(TARS) ELISA kit
Reactivity	Human
Range	12.5 pg/ml-800 pg/ml
Sensitivity	23.43 pg/ml
Tested Applications	ELISA
Target	TARS
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, cell lysates, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm