

Product Datasheet

Human HES1 ELISA Kit orb405433

Description	Human Transcription factor HES-1(HES1) ELISA kit
Reactivity	Human
Range	28 pg/ml-1800 pg/ml
Sensitivity	7 pg/ml
Tested Applications	ELISA
Target	HES1
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates, cell lysates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm