

## Product Datasheet

# Capsaicin Receptor (531-541) antibody (orb1129506)

<b>Description</b>	Rabbit polyclonal antibody to Capsaicin Receptor (531-541)...
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC-P
<b>Immunogen</b>	A synthetic peptide (YFSLKKEYVAS) of human capsaicin receptor protein (aa: 531-541) conjugated to KLH has been used as the immunogen.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1ug/ul
<b>Storage</b>	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
<b>Note</b>	For research use only.
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Filter purified
<b>Dilution Range</b>	Use at 1:1000 to 1:2000 dilution. This antibody has not been tested in other applications. The optimal dilution should be determined by the end user.
<b>Expiration Date</b>	12 months from date of receipt.