

Product Datasheet

Human TPP1 ELISA Kit orb405463

| | |
|----------------------------|--|
| Description | Human Tripeptidyl-peptidase 1(TPP1) ELISA kit |
| Reactivity | Human |
| Range | 12.5 pg/ml-800 pg/ml |
| Sensitivity | 3.12 pg/ml. |
| Tested Applications | ELISA |
| Target | TPP1 |
| Note | For research use only. |
| Expiration Date | Please enquire. |
| Application Notes | <ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates and cell lysates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm |