
Product Datasheet

Human TNFRSF21 ELISA Kit (orb405485)

Description	Human Tumor necrosis factor receptor superfamily member 21(TNFRSF21) ELISA kit
Reactivity	Human
Range	Request Information
Sensitivity	Request Information
Tested Applications	ELISA
Target	TNFRSF21
Note	For research use only.
Uniprot ID	O75509
Expiration Date	Please enquire.
Application Notes	• Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm