

Product Datasheet

Human TNFRSF25 ELISA Kit orb405486

Description	Human Tumor necrosis factor receptor superfamily member 25(TNFRSF25) ELISA kit
Reactivity	Human
Range	0.24 ng/ml-15 ng/ml
Sensitivity	0.06 ng/ml.
Tested Applications	ELISA
Target	TNFRSF25
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm