

Product Datasheet

Human UBB ELISA Kit orb405508

Description	Human ubiquitin B (UBB) ELISA kit
Reactivity	Human
Range	15.6 pg/ml-1000 pg/ml
Sensitivity	3.9 pg/ml.
Tested Applications	ELISA
Target	UBB
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm